

***Methods for the Analysis of Microcystins in
Water to sub-Parts Per Trillion Detection
Levels***

**Douglas Stevens,
Marijn Van Hulle, Gordon Fujimoto
Waters Corporation**

Introduction

- **Microcystins** are cyclic peptides produced by cyanobacteria
- They are **produced by overgrowth of algae**, especially at higher water temperatures – **Climate change, invasive species, global trade and agricultural practices can exacerbate** the problem
- Microcystins can be **toxic for plants and animals including humans**
- Once ingested, microcystins travel to the liver, via the bile acid transport system, where most is stored. Some remains in the blood stream and may contaminate tissue. Microcystins bind covalently to protein phosphatases thus disrupting cellular control processes. Their hepatotoxicity **may cause serious damage to the liver**. Microcystins can strongly inhibit protein phosphatases type 1 (PP1) and 2A (PP2A), and are linked to pansteatitis.
- Over 80 toxic variants are known

Introduction

- Microcystin-containing 'blooms' are a problem in countries worldwide including China, Brazil, Australia, the USA and much of Europe
- WHO action limit = **1000 ng/L (1ppb)** and methods are proposed including US EPA Methods 544 and 545 and European ISO 20179:2005(E) guideline
- For this reason sensitive detection is needed
- **This presentation will discuss the various approaches to sample introduction and detection suitable for proposed regulatory limits as well as the lower limits often desired for research purposes**

Experimental: UPLC conditions


- **Run Time: 7.50 min**
- **Inj Vol: 50.00 µL**
- Column: ACQUITY UPLC BEH C18 2.1x100mm, 1.7 µm
- Solvent A: 0.1% FA in 97/3 H₂O/ACN
- Solvent B: ACN
- Flow Rate: 400 µL/min
- Gradient:

Time (min)		%A	%B	Curve
0.00	100.0	0.0	Initial	
0.75	100.0	0.0	6	
5.00	20.0	80.0	6	
6.00	0.0	100.0	1	
7.50	100.0	0.0	1	

- Column Temp: 50.0 C


Experimental: MS Method Settings

- Time segmented MRM with Dwells to give 15 pts across each peak
- 2 MRMs per analyte
- Most microcystins singly charged except for DE-M-RR & M-RR



Performance Criteria

- The following performance criteria were evaluated
 - Linearity
 - Robustness
 - Sensitivity
 - Repeatability
 - Accuracy
 - Standard Addition



**Xevo TQ-S
shown with UPLC and APGC**

Mass Spectrometers


SQ, TQ, QTof, & QTof with Ion Mobility

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®




Results: Typical Chromatogram

- Bottom to top: Cylindrospermopsin (CYL), Anatoxin-A (ANA), M-DE-RR, M-RR, Nodularin (NOD), M-LR, M-YR, M-LY, M-LW and M-LF at 100 ng/L



Results: Sensitivity

- All analytes detected with good S:N at the 25 ng/L





Results: Linearity from 10 to 1000ng/L

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®


Compound name: Anatoxin-a
Correlation coefficient: $r = 0.999701$, $r^2 = 0.999402$
Calibration curve: $1978.2 \cdot x + 3045.67$
Response type: External Std, Area
Curve type: Linear, Origin: Include, Weighting: $1/x$, Axis trans: None

Anatoxin
 $R^2 = 0.9994$




Compound name: CYL
Correlation coefficient: $r = 0.998861$, $r^2 = 0.999722$
Calibration curve: $130.893 \cdot x + -183.659$
Response type: External Std, Area
Curve type: Linear, Origin: Include, Weighting: $1/x$, Axis trans: None

CYL
 $R^2 = 0.9997$





Compound name: Nodularin
Correlation coefficient: $r = 0.998806$, $r^2 = 0.999613$
Calibration curve: $70.7322 \cdot x + 6.27616$
Response type: External Std, Area
Curve type: Linear, Origin: Include, Weighting: $1/x$, Axis trans: None

Nodularin
 $R^2 = 0.9996$



Compound name: MC-YR
Correlation coefficient: $r = 0.999714$, $r^2 = 0.999429$
Calibration curve: $22.3124 \cdot x + -22.0103$
Response type: External Std, Area
Curve type: Linear, Origin: Include, Weighting: $1/x$, Axis trans: None

M-YR
 $R^2 = 0.9994$



Results: Robustness

- Due to lack of IS, signal stability was evaluated
- First and last (5th) 100 ng/L std were compared. Time difference between samples approximately 15.5 hrs. Areas within 10% for 5 of the 8 compounds, and 31% for the other 3

Compound	100 ng/L (1)	100 ng/L (5) +15.5h	ratio
M-YR	2061	2215	107%
M-LW	1361	1660	122%
M-LR	2770	3353	121%
M-LF	1522	1990	131%
NOD	6848	7380	108%
M-RR	28133	29911	106%
CYL	12027	12074	100%
ANA	190364	183937	97%

Results: Repeatability

- Short term repeatability evaluated using 5 replicates from the same vial at 100 ng/L
- % CV values for DW and surface water SW
- Repeatability is good with % CV values below 5%

Compound	% CV (n=5)	
	DW	SW
M-YR	4.6	3.0
M-LW	3.2	2.2
M-LR	2.2	1.8
M-LF	4.2	4.6
NOD	1.9	3.0
M-RR	1.6	0.5
CYL	5.0	2.7
ANA	1.6	0.5

Results: Accuracy

- Accuracy evaluated using spike/recovery experiments at 100, 200 and 400 ng/L in DW and SW
- DW recoveries are acceptable for compounds eluting in the middle of the run but too low for compounds eluting near the beginning and the end
- Suggests the need for IS to improve accuracy


% Spike Recovery Values

Name	DW 1	DW 2	DW 3	DW 4	DW 5
CYL	-54.7	-61.2	-59.6	-43.6	-40.8
Anatoxin-a	-2.4	-44.5	-47.5	-2.5	-59.3
MC-RR	7.9	5.0	7.4	3.7	3.0
Nodularin	4.0	3.5	5.9	3.8	6.6
MC-YR	-0.2	-1.8	2.5	2.5	9.4
MC-LR	2.7	-5.4	0.5	4.6	4.1
MC-LW	-20.6	-18	-29	-14.2	-19.6
MC-LF	-12.9	-14.9	-18	-13.1	-6.4

Results: Accuracy

- SW recoveries 30% - 150% and dependent on type of surface water
- Suggests the need for IS to improve accuracy
- **In the absence of an IS, standard addition can be used**
- Standard addition calculations can be done automatically in TargetLynx
- Following slides show quan of M-LR in Aquacheck Proficiency Test using double addition. Assigned value for M-LR, based on 14 results, is **720 ng/L**
- Sample analyzed using double standard addition:
 - Sample
 - Sample + 200 ng/L
 - Sample + 400 ng/L

Results: Standard Addition



Results: Standard Addition

- In the Sample List:
 - Use “Standard” for sample type
 - Add the spiked concentration
 - Assign a sample group for each set of samples

Spectrum Chromatogram Map Edit▼ Samples▼

	File Name	File Text	MS File	Inlet File	Bottle	Inject Volume	Sample Type	Conc A	Sample Group
1	1601926-01	AQ500-groep S32	test_JD	Microcysteines A1-B1	1:10	50.000	Standard		aquacheck
2	AQ_T200	AQ_T200	test_JD	Microcysteines A1-B1	1:11	50.000	Standard	200	aquacheck
3	AQ_T400	AQ_T400	test_JD	Microcysteines A1-B1	1:12	50.000	Standard	400	aquacheck

- In the TargetLynx processing method:
 - Check the option “Use Standard Addition?”

C:\MassLynx projects\20160204 Watergroep microcystines.PRO\20160204 stadand addition.qld - TargetLynx XS Method Editor

File Edit Update View Compound Help

Compound List


- 1: Anatoxin-a
- 2: Nodularin
- 3: CYL
- 4: MC-LF
- 5: MC-LR
- 6: MC-LW
- 7: MC-RR
- 8: MC-YR
- 9: MC-LY
- 10: DE-MC-RR

Calibration Properties

Compound Name	Value
Compound Name	MC-LR
Calibration Reference Compound	5: MC-LR
Concentration Units	
Concentration of Standard: Level	Conc A
Stock Concentration Factor	<input checked="" type="checkbox"/> 1.0000
Polynomial Type	Linear
Calibration Origin	Exclude
Weighting Function	1/X
Ignore Zero Level Standards?	<input checked="" type="checkbox"/> NO
Ignore Zero Level QCS?	<input checked="" type="checkbox"/> NO
Use Standard Addition?	<input checked="" type="checkbox"/> YES
Propagate Calibration Parameters?	<input checked="" type="checkbox"/> YES

Results: Standard Addition

- TargetLynx calculates endogenous concentration
- Reported conc. **724 ng/L in good agreement** with known value



Summary of Direct Injection High Sensitivity TQ Method

- Sensitive and robust method for 10 microcystins was developed:
 - 25 ng/L or lower LOD for all analytes
 - Signal stable for 15 hours during study
 - Repeatability <5% RSD
 - Linearity and %deviation excellent over the studied range
- Run time 7.5 min
- 50 µL of drinking water or surface water was injected directly
- In the absence of IS, standard addition yields accurate quantitation

Beyond PPB and PPT: PPQ Detection

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®

Google parts per quadrillion

All Images News Videos Shopping More Search tools

PPQ is an acronym for Parts per Quadrillion by allacronyms.com

Transports leading outdoor chemistries and particulate emissions of various indoor samples from human contact with...

PPQ means Parts per Quadrillion by allacronyms.com

Supplementary Material

FIG. 4

United States Patent Application Publication (18) Pub. No. US 2012/0325024 A1 (12) Pat. App. Pub. Date: Dec. 27, 2012 (51) Int. Cl. B01D 23/02 (2006.01) (52) U.S. Cl. 73/063.26 (57) ABSTRACT Impregnants are provided in the detection of atmospheric vapors by measuring their mass and their pressure, and analyzing them in situ. A novel detection limit of parts per quadrillion (ppq) concentrations are enabled by a combination of...

Sample Introduction

Carrier Gas

Sample Gas

Solid Surface with Residue

Ionization Source




Gas Flow

Reaction Region

Ion Movement

Mass Spectrometer

Parts per...
Parts per quadrillion (ppq)
- one gram of solute for every 10¹⁵ grams of solution (1 fg/g, 1 pg/kg).
- equivalent to:
• 1 gram of ink in Lake Ontario.
• 1 second in 31 million years
- Very few analytical techniques can measure with this degree of accuracy
- it is still used in some mathematical models of toxicology and epidemiology.







2D UPLC with 'At Column Dilution'


Loading Pump:

A pushes sample out of the autosampler

B dilutes the flow coming from the autosampler with aqueous to improve trapping



After the valve switch, flow from the column pump goes through the trap and onto the analytical column connected to the mass spec



Experimental Parameters

ACQUITY I-Class 2D-UPLC ACD Method

Trap Column
Oasis HLB Direct Conn
2.1 mm x 30 mm, 20

Analytical Column
Waters ACQUITY UPLC E
130Å, 1.7 µm, 2.1 mm X 50

UPLC System:
ACQUITY I-Class Bin
Solvent Manager (B
ACQUITY Column Manage

**Xevo
G2-XS
QToF**



250 µL Injections of
processed water sample

Purge Solution =
90/10 H₂O/MEOH

Wash Solution =
80/10 ACN/ IPA/MEOH/H₂O

Sample Temperature = 7 °C

Loading Pump

A = 67/33 Acetonitrile/Water

B = Water

A & B have 0.3% Formic Acid

Gradient (Flow=0.80 mL/min)

<u>Time</u>	<u>A%</u>	<u>B%</u>		<u>Time</u>	<u>A%</u>	<u>B%</u>
0.0 min	15	85	Switch at 2.25	0.0 min	90	10
3.3 min	15	85		3.3 min	90	10
3.4 min	95	5				
6.0 min	95	5				
6.5 min	15	85				
				7.3 min	5	95
				8.5 min	5	95
				9.0 min	90	10
			Switch at 10.25			
11.0 min	15	85		11.5 min	90	10
11.5 min	15	85				

Column Gradient Pump

A = Water with

B = Acetonitrile with

A & B have 0.01% Formic Acid

Gradient (Flow=0.45 mL/min)

t=0: Flow from autosampler goes through trap column to waste
t=2.25: valve switch to send flow from the column pump through the trap and onto the analytical column
t=10.25 switch back to t=0 flow path

Loading Pump:

A is used to push the sample out of the autosampler


B is used to dilute the flow coming from the autosampler with aqueous before it gets to the trap column

Blank Injections 1 - 4

Water Extraction Blank 1 Acquired 16-Dec-2015 / 02:44:14

Test_MCLR_601 Sm (Mn, 1x2)


1: TOF MS ES+
995.557 0.0200Da
1.06e3



Water Extraction Blank 3 Acquired 16-Dec-2015 / 03:12:23

Test_MCLR_603 Sm (Mn, 1x2)


1: TOF MS ES+
995.557 0.0200Da
1.36e3



Water Extraction Blank 2 Acquired 16-Dec-2015 / 02:58:18

Test_MCLR_602 Sm (Mn, 1x2)


1: TOF MS ES+
995.557 0.0200Da
1.62e3



Water Extraction Blank 4 Acquired 16-Dec-2015 / 03:26:27

Test_MCLR_604 Sm (Mn, 1x2)

1: TOF MS ES+
995.557 0.0200Da
1.53e3



Blank Injection 5 and & RO Water Injection 1


Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Water Extraction Blank 5 Acquired 16-Dec-2015 / 03:40:32

Test_MCLR_605 Sm (Mn, 1x2)

1: TOF MS ES+
995.557 0.0200Da
1.60e3




RO 1 ppt Extraction 1

Acquired 16-Dec-2015 / 03:55:02


Test_MCLR_606 Sm (Mn, 1x2)

1: TOF MS ES+
995.557 0.0200Da
1.34e4
Area



Test_MCLR_606 Sm (Mn, 1x2)

S/N:RMS=64.71 1: TOF MS ES+
995.557 0.0200Da
1.34e4




RO Water Injections 2 & 3

RO 1 ppt Extraction 2 Acquired 16-Dec-2015 / 04:09:05


Test_MCLR_607 Sm (Mn, 1x2)

1: TOF MS ES+
995.557 0.0200Da
1.80e4
Area



Test_MCLR_607 Sm (Mn, 1x2)


S/N:RMS=90.10 1: TOF MS ES+
995.557 0.0200Da
1.80e4



RO 1 ppt Extraction 3 Acquired 16-Dec-2015 / 04:37:37


Test_MCLR_608 Sm (Mn, 1x2)

1: TOF MS ES+
995.557 0.0200Da
1.79e4
Area



Test_MCLR_608 Sm (Mn, 1x2)

S/N:RMS=75.72 1: TOF MS ES+
995.557 0.0200Da
1.79e4




RO Water Injections 3 & 4

RO 1 ppt Extraction 4 Acquired 16-Dec-2015 / 04:51:40

Test_MCLR_609 Sm (Mn, 1x2)


1: TOF MS ES+
995.557 0.0200Da
1.81e4
Area



Test_MCLR_609 Sm (Mn, 1x2)

S/N:RMS=73.83


1: TOF MS ES+
995.557 0.0200Da
1.81e4



RO 1 ppt Extraction 5 Acquired 16-Dec-2015 / 05:05:43

Test_MCLR_610 Sm (Mn, 1x2)


1: TOF MS ES+
995.557 0.0200Da
1.74e4
Area



Test_MCLR_610 Sm (Mn, 1x2)

S/N:RMS=67.51

1: TOF MS ES+
995.557 0.0200Da
1.74e4



Instrument Detection Limits

Five Analyses of RO Water Samples Spiked at 1 ppt

<u>Injection</u>	<u>RT</u>	<u>Area</u>	<u>RMS S:N</u>
1	5.48	355.2	64.71
2	5.48	527.7	90.10
3	5.48	525.8	75.72
4	5.48	500.8	73.83
5	5.48	527.8	67.51

Average =	487.5	74.37
Std Dev =	74.8	9.87
RSD =	15.3%	13.3%


Translates to LOD of approximately <100ppq

Spectrum for 1 ppt of Microcystin LR

10 mL Water MCLR Spike

Test_MCLR_562 1007 (5.490) Cm (1003:1007-989:993)

1: TOF MS ES+
1.06e5




Why Use a QToF?

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®

Item name: Water_005
Channel name: 1: TOF MSe (50-1200) 6eV ESI+

Because you need to know more.



Linearity for Microcystin LR

Compound name: Microcystin LR


Correlation coefficient: $r = 0.999008$, $r^2 = 0.998017$

Calibration curve: $197.411 * x + -5.88617$

Response type: External Std, Area

Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None

Calibration Curve for Microcystin LR 0.010ppb to 100ppb



Summary of 2D LC QToF Method

- Sensitive (PPQ) detection for MC-LR was achieved using 2D LC
 - 100 pg/L detection limit achieved
 - Linearity from 0.01 to 100ppb $>0.99 R^2$
 - Matrix effects were evaluated – minimal and acceptable
- Run time <15 min
- **Overall Summary**
 - Fit for purpose LOD can be achieved with high performance TQ through direct injection and little or no sample pretreatment
 - 2D LC and SPE provide enrichment factors leading to fit for purpose LOD for lower sensitivity TQ or unparalleled sensitivity for research purposes when paired with high performance MS

Acknowledgements

- De Watergroep, Belgium



- Ministry of Environment and Climate Change, Toronto

