

Basics and Importance of Total Organic Carbon (TOC) Measurements in Environmental Samples



Dan Davis

Shimadzu Scientific Instruments

Product Manager for TOC/AA/ICP/X-Ray


TOC is the measure of organic carbon in a sample

- **measures all carbon converted to CO_2**
- **not source specific**
- **reacts with disinfectants to form THMs**




Organic carbon is derived from living molecules

- **Sugars**
- **Metabolic Activities**
- **Decaying Vegetation**
- **Carbohydrates etc.**



Inorganic Carbon is non-living carbon

- $\text{CO}_{(2)}$ [dissolved]
- Carbonates
- Bicarbonates
- Carbonyls etc.



Carbon Glossary

TC		
IC	TOC	
	POC*	NPOC
		DOC
		PtOC

- TC = Total Carbon
- IC = Inorganic Carbon
- TOC = Total Organic Carbon
- POC = Purgable Organic Carbon
- NPOC = Non-Purgable Organic Carbon
- DOC = Dissolved Organic Carbon
- PtOC = Suspended (Particulate) Organic Carbon



Necessary definitions for TOC analyses

Purgable Organic Carbon

- organic carbon removed from an acidified sample by purging
- referred to as Volatile Organic Compounds (VOC)

Non-Purgable Organic Carbon

- organic carbon remaining in an acidified sample after purging

Dissolved Organic Carbon

- organic carbon remaining in a sample after filtering the sample, typically using a 0.45 micrometer filter


Suspended Organic Carbon

- the carbon in particulates too large to pass through a filter



Estimating Organic Carbon

Various methods are used to estimate organic carbon in environmental samples



BOD



**Biological Oxygen Demand
measures oxygen consumption by
bacteria over 5 days**

**Also measures oxidation of NH_3 to
 NO_2 and NO_3 ,**



COD

- **Chemical Oxygen Demand measures the oxygen required for chemical oxidation**
- **Organic matter, sulfide, iron, nitrogen all oxidize**
- **Some organic molecules do not oxidize**



COD and TOC

There is a theoretical correlation between COD and TOC

$$\text{Theoretical COD} = \text{TOC} \times \text{O}_2/\text{C}$$

$$\text{COD} = \text{TOC} \times 32/12$$

$$\text{COD} = \text{TOC} \times 2.67$$



TOC and DBPs

Minimizing Disinfection byproducts requires TOC removal

Source Water TOC (mg/L)	Source Water Alkalinity (mg/L as CaCO ₃)		
	0 – 60	60 – 120	Greater than 120
2.0 – 4.0	35.0 %	25.0 %	15.0 %
4.0 – 8.0	45.0 %	35.0 %	25.0 %
Greater than 8.0	50.0 %	40.0 %	30.0 %



Some CWA permits require TOC



TOC and Standard Methods

Technique	Method
High Temperature Oxidation	SM 5310B
Wet Chemical Oxidation	SM 5310C



To measure TOC you must have an Analyzer

- Measure organic and/or inorganic carbon in water or soil
- Used for DBP rule to estimate THM formation
- Correlate TOC with BOD and COD



Measuring TOC



Most TOC analyzers remove TIC and measure NPOC

- NPOC = remove TIC with acid and measure remaining Carbon
- $TOC = TC - TIC$
 - Measure TC and TIC and subtract
 - Prone to error
 - Preferred for samples that foam




Detecting CO₂

- IR detection
 - requires carrier gas
 - specific to CO₂
- Membrane conductivity detection
 - Works by gas diffusion
 - Anything that passes through is detected
 - No carrier gas required



NDIR detection is practically interference free



Maximum Absorbance
2340 cm⁻¹



TOC Oxidation Method

TOC analyzers oxidize carbon by high temperature or by chemical oxidation

- High Temperature
 - Combustion
 - Catalyst assisted oxidation
- Chemical Oxidation
 - Persulfate
 - Catalyzed by heat
 - UV irradiation
 - Heat and UV irradiation



TOC is an Oxidation Reaction


**When thinking about oxidation of TOC
remember that:**

**For a reaction to happen a collision must
occur...**



TOC is an Oxidation Reaction

Combustion is a rapid form of oxidation that “burns” compounds to CO_2 and H_2O



TOC is an Oxidation Reaction

For the reaction, we know:

- Oxygen and “fuel” must collide (vapor state)
- More oxygen = greater chance to collide
- Higher temperature = more collisions
- More time = more collisions



Combustion TOC Analyzer methods

Combustion TOC Analyzer methods are carried out in combustion tubes at 950°C or above

- Carrier gas is pure oxygen
- Detector is usually IR but can be Coulometric
- Combustion of carbonates...soil, rocks, or coal



Catalytic Combustion




Catalytic Combustion oxidizes organic compounds to CO_2 and H_2O at lower temperatures



- Samples are injected onto a metal, or metal oxide catalyst where collisions occur
- Reaction steps involve adsorption and desorption
- Organic compounds react with molecular oxygen on the surface of the catalyst
- Oxidation efficiency is a function of catalyst surface area and catalyst contact time



TOC is an Oxidation Reaction



TOC is an Oxidation Reaction



TOC is an Oxidation Reaction

- Oxidation of larger molecules proceeds stepwise as the organics travel through the catalyst.
- Oxidation efficiency \sim temperature (surface area / flow rate).
- As the surface area, flow rate, temperature increases so does the chance of collisions.
 - as temperature rises unwanted side reactions also increase



Why Temperature Matters

- Melting point of NaCl = 801°C
- $\text{NaCl} + \text{H}_2\text{O} \rightarrow \text{Na}_2\text{O} + 2\text{HCl}$
- $\text{Na}_2\text{O} + \text{SiO}_2 \rightarrow \text{Na}_2\text{SiO}_3$



Why Temperature Matters

Catalytic Combustion TOC methods are carried out in combustion tubes at $\geq 680^{\circ}\text{C}$

- Lower temperature increases detector life
- Lower temperature increases combustion tube life
- Catalytic Combustion is often used to determine TOC in drinking water, wastewater, and seawater.



Chemical Oxidation

Chemical oxidation reacts an oxidizing chemical with organic compounds to CO_2 and H_2O

- Homogenous reaction usually acid or basic persulfate
- Catalyst speeds the reaction
 - Heat
 - UV
- The catalyst converts $\text{S}_2\text{O}_8^{2-}$ to the SO_4^\bullet radical
 - Stable persulfate + heat = free radical
 - The highest demand for radical is the matrix
- The number of collisions is proportional to temperature
 - Persulfate decomposes in high heat



Oxidation Techniques: Advantages and Disadvantages

Technique	Advantages	Disadvantages
Combustion (HTC)	High efficiency Soil, coal, solids	High detection limits Hazardous temperatures
Catalytic Combustion (HTCC)	High Efficiency Lower temperature than HTC Tolerates High Salt	Higher detection limit than PO Memory effects of catalyst Catalyst poisoning
Chemical Oxidation (PO)	Low Detection Limits No memory effects	Chloride interferes Lower oxidation efficiency



Auto Samplers Considerations



- High Sample Capacity
- Septum Piercing
- Stirring

- Lower Sample Capacity
- Fits Any Size Vial
- Stirring




Total Nitrogen Option



- TN add-on for TOC System
- Simultaneous Measurements of Carbon & Nitrogen
- 5 minute analysis



Options for Solid Sampling



- TC & IC Measurements
- Pt/CoO₃ Catalyst
- Adjustable Temperature



Online Options



- On-Line Continuous Monitoring
Carbon, Nitrogen, Phosphorus
- Multi-Steam Particulate Sampler
Up to 6 streams
- Remote Start, Stop and Alarms
- Same method as lab TOC



Northeast Environmental Seminar

November 6th

Shimadzu Scientific Instruments Headquarters, Columbia MD

Information at <http://www.ssi.shimadzu.com/news/>

Speakers from EPA, Restek, EST, Inorganic Ventures



Thankyou!

Any Questions?



www.ssi.shimadzu.com